

Biax Experiment

Exp. Name: p5483WGFracNSPPosc
Operator(s): Wood, Borate

Date/Time: 30 Jan. 2021
Hydraulics start: 3808.1
Hydraulics end: 3846.5

Sample Block Thickness w/ no gouge:

Layer Thickness (total on bench): mm @sample

Under Load: mm

Material (Qtz, Granite, ?): WG, Ex-situ Fractured (sample p5346).

Particle Size, Size Distribution :

Load Cells:

Contact Area: 0.0022292545 m²

Load cell name	Calibrations (mV/kN)	Target stress (MPa)	Init. Voltage	Volt. @ load
44mm Solid Horiz	119.3033 (V/MPa): 0.26596	2, 9.25, 11, 13, 15, 18	0.999	1.53091, 3.45911, 3.92453, 4.45645, 4.98836, 5.78623

Vessel Pressures:

Pore Fluid: DI H2O

	Calibrations (V/MPa)	Pressures (MPa)	Init. Voltage	Volt. @ load
Pc	Gain: 0.1456	8.25, 10.5, 12, 13.5, 12	-0.1495	1.05137, 1.37888, 1.59722, 1.81557, 1.59722
PpA	1.5177	2.6	-0.114	3.83197
PpB	1.483	2.6, 1.4	-0.56	3.29585, 1.51623

Data Logger Used: 16 channel

Horiz. DCDT: short rod

0.6545 mm/V

Control File: No

Vert. DCDT: Trans-Tek 2

2.8498 mm/V

Purpose/Description: p5369 experiment follow-up.

DAET: oscillate PP and limited NS oscillations and increasing effNS (modulate fracture closure).

Active-source ultrasonic recording: 8 transmitters, 10 receivers using SDS sideblocks v2.

Acoustics Blocks used: SDS L-block v2

Horiz. Servo Settings				Vert. Servo Settings	
P	500	Datten	10	P	Datten
I	80	Feedback	512	I	Feedback
D	10	E-gain	700	D	E-gain

@ Hyd. Power Supply (HPS)	Chilled water at HPS	Chiller Unit	Process water at Chiller
14. Tank Temp. (C)	1. Temp. In (F)	6. Panel Temp. (F)	10. Temp. In (F)
15. Temp. Out (C)	2. Pres. In (psi)	7. Panel Pres. (psi)	11. Pres. In (psi)
16. Pres. Out (psi)	3. Temp. Out (F)	8. Near Pres. In (psi)	12. Temp. Out (F)
	4. Pres. Out (psi)	9. Near Pres. Out (psi)	13. Pres. Out (psi)
	5. Flow (lpm)		

Experiment Notes

```
# Ppa, Ppb disp and load are switched in mech file
# 80 Int disp offset
# 5600 effNS @ 10MPa
# 15577 ↗ 100Hz, Ppa,Ppb ↗ 2.6 MPa. Ppb ↗ 1.2MPa. Run1, Run2, Run3
# 504550 effNS @ 12.5 MPa
# 504900 ↗ 100Hz, Ppa,Ppb ↗ 2.6 MPa. Ppb ↗ 1.2MPa. Run4, Run5, Run6
# 859600 effNS @ 15 MPa
# 859900 ↗ 100Hz, Ppa,Ppb ↗ 2.6 MPa. Ppb ↗ 1.2MPa. Run7, Run8
# 1140000 ↘ 1Hz, refill Ppa, empty Ppb, ↗ 100Hz, Run9
# 1336600 ↘ 1Hz, Ppa & Ppb ↘ 0MPa, Pc ↘ 0MPa, NS ↘ 1MPa
# 1398660 effNS @ 17.5 MPa
# 1399880 ↗ 100Hz, Ppa,Ppb ↗ 2.6 MPa. Ppb ↗ 1.2MPa. Run10, Run11, Run12
# 1778000 ↘ 1Hz
# 1783965 effNS @ 20 MPa
# 1784250 ↗ 100Hz, Ppa,Ppb ↗ 2.6 MPa. Ppb ↗ 1.2MPa. Run13, Run14, Run15
# Labview crashed... memory overflow...begin Part2
```

```
# 300 ↗ 100Hz, Ppa,Ppb ↗ 2.6 MPa. Ppb ↗ 1.2MPa. Run16 - 19
# 393000 ↘ 1Hz, effNS @ 12.5MPa
# 397000 ↗ 100Hz, Ppa,Ppb ↗ 2.6 MPa. Ppb ↗ 1.2MPa. Run20 - 23
# 640633 redo Run22 after 0.4MPa amp oscillation (2nd in set of 7)
# 655000 verasonics computer turned off, maybe just screen
# 640633 redo Run22 after 1Hz
```